



dental spa

user manual



visit us at www.violight.com



Violight Dental Spa, with cleaning cup and contour tray

see our other UV products at www.violight.com



Our award-winning countertop and personal toothbrush sanitizers.



1-year limited warranty

Violight Inc. warrants to the original purchaser of this product that it is free from defects in materials and workmanship for a period of one year from date of purchase. If the unit becomes defective within the one-year period, call 1-800-950-6101 with original purchase receipt.

Within one year of purchase, VIOlight Inc. will repair or replace the product we deem defective at no charge, provided the product has not been abused, misused, or damaged, and was used according to these instructions.

This warranty gives you specific legal rights and you may also have other rights, which vary from state to state.

For return address information and instructions, go to www.violight.com or call 1-800-950-6101 or email customerservice@violight.com.

copyright ©2007 Violight Inc.



dental spa

Congratulations on your purchase of the VIOlight Dental Spa. As a health-conscious consumer, you're aware of the benefits of keeping your dental products clean and sanitized.

With your new VIOlight Dental Spa, keeping your dentures, orthodontics, mouthguards, and whitening trays is easier than ever. With the push of a button, your Dental Spa does the work for you! See the blue glow, feel the sonic vibration, and you'll know it's working! Within 5 minutes, the germicidal ultraviolet light will kill up to 99.9% of the germs it comes into contact with, and the sonic action will remove particles and residue for sparkling clean results.

By following these instructions and safeguards, you and your family will enjoy your Dental Spa for many years to come.



dental spa

using your dental spa:

Before use, unscrew the battery door compartment and insert 3 "AA" batteries.

1. Push the small white clip to open the lid of your Dental Spa. Remove the cup and contour tray from the inside compartment by lifting straight up and out.
2. Set the curved tray aside.
3. Fill the cup with water to the indicator line.
4. Place the cup back into the Dental Spa. When seated, push down firmly for a tight fit.
5. Place items to be cleaned directly into the cup so they are submersed in the water. If you are cleaning smaller items, you can place them on the curved tray and then insert into the cup of water.
6. Close the lid of the Dental Spa.
7. **For UV + Sonic Cleaning:** Press the [UV + Sonic] button and a 5-minute cycle will start. The LED indicator light will flash red and blue. The Dental Spa will turn off automatically.

For UV Cleaning Only: Press the [UV Only] button and a 5-minute cycle will start. The LED indicator light will flash blue. The Dental Spa will turn off automatically.
8. **Manual Shut-Off:** To interrupt either cycle, simply press the appropriate button and hold for 5 seconds.
9. After your Dental Spa cycle is complete, remove your items and discard the water from the cup. Rinse and dry the cup so it's ready for the next use!

safeguards:

- Children should be properly supervised.
- Do not place unit where it can be in danger of falling into sinks, toilets or bathtubs.
- Always make sure your Dental Spa lid is closed before using.
- Do not look directly at the lit bulb.
- If the unit is dropped or damaged, or if the unit remains on beyond the 5-minute sanitizing cycle, please call or email customer service.

cleaning & care:

- Remove and rinse the internal cup and contour tray after each use.
- NEVER rinse the base, place it under running water, or immerse it for any reason. To clean, simply wipe with a moist cloth or sponge as needed.

if there's a problem:

check the following:

- Did you press the blue buttons hard enough?
- Are the batteries in need of replacement?

If the problem persists, please call VIOlight directly — not the catalog or store from which you purchased. We know your VIOlight unit better than anyone, and we're here to help you: send email to customerservice@violight.com or call 888-882-2367. Visit us on the web at www.violight.com.